

PEDIATRIC SURGERY

Circumcision

Preference Card

1. Put in dorsal penile nerve block and ring block with marcaine.
2. Retract foreskin and separate from glans using a lacrimal probe.
3. Divide frenular artery using a mosquito and bovie.
4. Size the Plastibell then cut a dorsal slit by clamping with a straight clamp, releasing then cutting with bovie.
5. Place the Plastibell on the glans within the foreskin and tie in place.
6. Cut just above the tie with bovie then remove Plastibell. Coagulate any visible vessels in the sub-Q of the penis.
7. Sew cut edges of foreskin together with 5-0 plains. Place at 12, 3, 6, 9 o'clock and tag then place 3 between each of these 4.
8. Dress with neosporin, xeroform and 2x2's covered by a diaper.